Uranyl Acetate Staining protocol in the IMF

Uranyl acetate is in a 2% solution.

- 1. Use the designated pipettes and Radiation shield stored in the locked drawer near the fumehood in IMF SE2-90G
- 2. Lab coat, gloves, and safety glasses are required when handling uranyl acetate
- 3. Prepare samples in the workbench designated for radioactive materials located in IMF SE2-90G
- 4. Remove TEM grid from its grid storage using the TEM specified tweezers.
- 5. Place TEM grid in your grid case
- 6. Apply 4 μl of sample to the grid and wait 5-10 min
- 7. Cut out a 3x3 inch piece of parafilm and place within the specified radioactive area. place 3 drops (40 μ l each) of HPLC water on the top of the parafilm.
- 8. Remove the uranyl acetate from the locked drawer located in IMF room SE2-90G and place in the specified radioactive area.
- 9. Place 3 drops (40 ul each) of uranyl acetate in the lower part of the parafilm.
- 10. Remove TEM grid from case and wash in HPLC water (10 seconds each droplet). Then wash in uranyl acetate (10 seconds each droplet).
- 11. Store back in grid case and let sit for 5-10 min.
- 12. Dry using a kimwipe by placing the tip of the kimwipe next to the grid.
- 13. Let dry overnight
- 14. Place uranyl acetate back in drawer and make sure the drawer is locked.

Alternative: TEM protein sample staining with uranyl acetate (ver.1, 2023.06.12)

We use carbon coated formvar-copper grids (CF400-CU-50, Electron Microscopy Sciences).

- 1. The grids are treated with a plasma glow discharger to make the grid surface hydrophilic.
- Sample solution (typically use 5 μL, 0.05 0.1 mg/mL protein) is applied to the carbon coated side of the treated grid and incubated for 30 seconds. Excess liquid was blotted with filter paper. Wait about 5 min to make the surface dry.
- 3. To remove salt from the grid surface, ultra-pure water (5 μL) is applied to the grid surface and water is immediately wicked away with filter paper. Wait about 5 min to make the surface dry. If necessary, repeat this process one more time.
- 4. Grid is then stained with 5 μ L 1 2 % uranyl acetate and excess stain was wicked away with filter paper. Wait about 5 min to make the surface dry. The sample grid is ready for TEM observation